



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A34730**

PRODUCT:	Troponin I Monoclonal Antibody
SOURCE:	Mouse ascites fluid
IMMUNOGEN:	Purified human cardiac Troponin I subunit.
FORM:	Liquid
PURITY:	>95% by SDS-PAGE
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% NaN <sub>3</sub> .
SUBCLASS:	IgG
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.