



PRODUCT SPECIFICATION

CATALOG NUMBER: A34750

PRODUCT:	Troponin C Monoclonal Antibody
IMMUNOGEN:	Human cardiac Troponin C
FORM:	Purified from murine ascites
PURIFICATION:	Purified by chromatography on protein A sepharose. Purity is tested by electrophoresis.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% of NaN ₃
ISOTYPE:	IgG _{2a}
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.