



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A34760**

PRODUCT:	Troponin C Monoclonal Antibody
IMMUNOGEN:	Human cardiac Troponin C
FORM:	Purified immunoglobulin from murine ascites
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity is tested by electrophoresis.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% NaN <sub>3</sub>
ISOTYPE:	IgG1
STORAGE CONDITIONS:	2-8°C