



PRODUCT SPECIFICATION

catalog number:	A34770
product description:	Troponin C (TnC) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human cardiac troponin C
form:	Liquid
purification method:	Protein A chromatography
buffer:	PBS, pH 7.4
preservative:	0.09% NaN ₃
isotype:	IgG _{2b}
concentration:	[lot specific] mg/mL
specificity:	This antibody reacts with human TnC.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.