



## PRODUCT SPECIFICATION

catalog number:	<b>A34780</b>
product description:	Troponin I (137-148) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Bioreactor cell culture supernatant
form:	Frozen liquid
purification method:	Protein G chromatography
purity:	>98% by HPLC
buffer:	PBS, pH 7.3
preservative:	none
isotype:	IgG <sub>1</sub> , kappa
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	-20°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.