



PRODUCT SPECIFICATION

CATALOG NUMBER: A34790

PRODUCT:	Troponin I Monoclonal Antibody
IMMUNOGEN:	Purified human cardiac Troponin I
FORM:	Liquid
PURIFICATION:	Immunoglobulin purified by protein G sepharose.
BUFFER:	PBS, pH 7.3
PRESERVATIVE:	None
ISOTYPE:	IgG ₁ /K
STORAGE CONDITIONS:	2-8°C short term, -20° C long term
APPLICATION:	ELISA, Western Blot
COUNTRY OF ORIGIN:	Canada
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.