



PRODUCT SPECIFICATION

CATALOG NUMBER: A34800

PRODUCT:	Troponin I Monoclonal Antibody
SOURCE:	Murine Hybridoma
IMMUNOGEN:	Purified human cardiac Troponin I
FORM:	Liquid
PURIFICATION METHOD:	Protein G
BUFFER:	PBS, pH 7.3
PRESERVATIVE:	0.05% Sodium Azide
ISOTYPE:	IgG _{1/K}
STORAGE CONDITIONS:	2-8 ⁰ C short term -20 ⁰ C long term
APPLICATION:	ELISA or Western Blot
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.