



## PRODUCT SPECIFICATION

catalog number:	<b>A34803</b>
product description:	TNF Alpha Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Recombinant TNF alpha
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	BBS 0.1M pH 8.0
preservative:	None
isotype:	IgG1
concentration:	[lot specific]
specificity:	Recognizes native and recombinant Human TNF alpha
pairing:	This antibody can be used in sandwich ELISA as the capture antibody
storage:	2-8°C
expiration date:	[lot specific]
hazard/biohazard:	There is no known biohazard associated with this product.