



PRODUCT SPECIFICATION

catalog number:	A34804
product description:	TNF Alpha Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Recombinant TNF alpha
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	PBS 0.1M, pH7.4
preservative:	0.5g/L NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific]
specificity:	Recognizes native and recombinant Human TNF alpha
pairing:	This antibody can be used in sandwich ELISA as the capture antibody
storage:	2-8°C
expiration date:	[lot specific]
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.