



PRODUCT SPECIFICATION

CATALOG NUMBER: A34810

PRODUCT:	Troponin I Monoclonal Antibody
IMMUNOGEN:	Purified human cardiac Troponin I subunit.
FORM:	Liquid
PURITY:	≥ 95% by SDS-PAGE.
BUFFER:	0.1M Tris-Glycine, pH 8.0
PRESERVATIVE:	0.1% NaN ₃
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Human cardiac specific Troponin I, a.a. 59-79.