



## PRODUCT SPECIFICATION

catalog number:	<b>A34910</b>
product description:	Skeletal Troponin I Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified human fast skeletal troponin I
source:	Bioreactor cell culture supernatant
form:	Frozen liquid
purification method:	Protein G chromatography
buffer:	PBS, pH 7.3
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub> , kappa
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.