



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A34920**

ANTIBODY TO:	Skeletal Troponin I
SOURCE:	Murine Hybridoma
IMMUNOGEN:	Purified human fast skeletal troponin I
FORM:	Liquid
PURIFICATION METHOD:	Protein G
BUFFER:	PBS, pH 7.3
ISOTYPE:	IgG <sub>1</sub> /K
STORAGE CONDITIONS:	-20°C
APPLICATION:	ELISA, Western Blot
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.