



PRODUCT SPECIFICATION

catalog number:	A34992
product description:	FABP Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
immunogen:	Purified Human FABP3
form:	Liquid
purification method:	Protein A Affinity Purified
purity:	≥ 95%
buffer:	50mM NaCl, 20 mM PBS pH7.4
preservative:	0.05% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific]
specificity:	Specific for Human FABP in Double Antibody Test
expiration date:	[lot specific]
storage:	2-8°C short term. For long term this can be stored at -20°C
pairing:	This antibody forms a suitable pair with clone A34991
hazard/biohazard:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known hazard or biohazard associated with this product