



## PRODUCT SPECIFICATION

catalog number:	<b>A34995</b>
product description:	FABP Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified human native fatty acid binding protein
form:	Liquid
purification method:	Protein A affinity chromatography; 0.2µm Filtered
purity:	>98% by SDS-PAGE
buffer:	PBS, pH 7.2
preservative:	none
concentration:	[lot specific] mg/mL
storage:	-20°C
expiration:	2 years from date of receipt
hazard/biohazard:	There is no known biohazard associated with this product.