



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A35230**

PRODUCT:	Human Chorionic Gonadotropin, Intact
IMMUNOGEN:	hCG from human urine
FORM:	Chromatographically purified immunoglobulin fraction in 0.05 M PBS, pH 7.4
PRESERVATIVE:	0.1% NaN <sub>3</sub>
AFFINITY CONSTANT:	1 x 10 <sup>10</sup>
SUBCLASS:	IgG
STORAGE CONDITIONS:	2-8°C