



PRODUCT SPECIFICATION

CATALOG NUMBER: A36020

PRODUCT:	Human Growth Hormone (HGH) Monoclonal Antibody
IMMUNOGEN:	Purified human growth hormone
BUFFER:	20mM sodium phosphate, 0.15M NaCl, pH 7.2
PURIFICATION:	Protein A
PURITY:	>95%
PRESERVATIVE:	0.05% NaN ₃
SUBCLASS:	IgG _{1,k}
STORAGE CONDITIONS:	2-8°C short term, -20°C long term