



PRODUCT SPECIFICATION

catalog number:	A36025
product description:	hGH Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Pituitary GH
source:	Murine ascites
form:	Liquid; Purified IgG
purity:	>95% by SDS-Gel
buffer:	PBS, pH 7.0
preservative:	0.02% NaN ₃
isotype:	IgG _{2b}
concentration:	[lot specific] (A280nm, 1mg/ml, 1cm=1.4)
storage:	2-8°C short term; -20°C long term
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.