



PRODUCT SPECIFICATION

CATALOG NUMBER: A36040

PRODUCT:	Human Growth Hormone (HGH) Monoclonal Antibody
IMMUNOGEN:	Purified human growth hormone
PURIFICATION METHOD:	Protein A
PURITY:	>90% by SDS-PAGE
BUFFER:	10 mM phosphate with 150 mM NaCl, pH 7.4
PRESERVATIVE:	0.1% Sodium Azide
SUBCLASS:	IgG _{1,k}
STORAGE CONDITIONS:	2-8°C short term, -20°C long term Avoid repeated freeze-thaw cycles.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.