



PRODUCT SPECIFICATION

CATALOG NUMBER: A36050

PRODUCT:	Human Growth Hormone (HGH) Monoclonal Antibody
IMMUNOGEN:	Purified human growth hormone
FORM:	Purified immunoglobulin from Mouse Ascites
PURIFICATION METHOD:	Protein A
PURITY:	>90% by SDS-PAGE
BUFFER:	0.01M phosphate, 0.150 M NaCl, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG _{1,Kappa}
STORAGE CONDITIONS:	2-8°C short term, -20°C long term Avoid repeated freeze-thaw cycles.
PART OF A COOPERATIVE PAIR:	Yes, A36030
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.