



PRODUCT SPECIFICATION

CATALOG NUMBER: A36130

PRODUCT:	<u>Staphylococcus aureus</u> Enterotoxin Monoclonal Antibody
SOURCE:	Immunoglobulin from mouse ascites
IMMUNOGEN:	<u>Staphylococcus aureus</u> Enterotoxin A
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity is tested by electrophoresis.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG _{2a}
STORAGE CONDITIONS:	4°C
APPLICATION:	Detection of staphylococcal enterotoxins in ELISA.
REACTIVITY:	Staphylococcal enterotoxins B, C1, C2 and D