



## PRODUCT SPECIFICATION

catalog number:	<b>A36330</b>
product description:	HbA1c Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Glycated (in vitro) synthetic peptide corresponding to the N-terminal of beta-globin
source:	Murine
form:	Liquid
purification method:	Protein G chromatography
purity:	>95% by SDS-PAGE
buffer:	PBS, pH 7.4 ±0.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
cross-reactivity:	This clone specifically recognizes HbA1c and does not cross react with HbA0.
pairing:	This antibody forms a match pair with Hemoglobin antibody A36340
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.