



## PRODUCT SPECIFICATION

catalog number:	<b>A36340</b>
product description:	Hemoglobin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified human hemoglobin
form:	Liquid
purification method:	Protein G chromatography
purity:	>95% by SDS-PAGE
buffer:	PBS, pH 7.4 ±0.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
cross-reactivity:	This antibody cross-reacts with human. It does not cross-react with bovine, chicken, horse, porcine or sheep.
pairing:	This antibody forms a match pair with HbA1c antibody A36330.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.