



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A36440**

PRODUCT:	Human Hemoglobin Monoclonal Antibody
IMMUNOGEN:	Human Hemoglobin
PURITY:	>95% by SDS-Page
FORM:	Purified IgG from mouse ascites
BUFFER:	Phosphate buffer saline, pH 7.0
PRESERVATIVE:	0.02% NaN <sub>3</sub>
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8°C short term; -20°C long term
PAIRING INFORMATION:	This antibody form a suitable sandwich when paired with A36430.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.