



PRODUCT SPECIFICATION

CATALOG NUMBER: A37010

PRODUCT:	Human IgE Monoclonal Antibody
IMMUNOGEN:	Human IgE (serum)
FORM:	Chromatographically purified immunoglobulin from mouse ascites.
BUFFER:	0.05M PBS, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
AFFINITY CONSTANT:	3 x 10 ¹⁰
SUBCLASS:	IgG _{2b}
STORAGE CONDITIONS:	2-8°C (Do not freeze)
CROSS REACTIVITY:	IgG < 1%, IgA < 1%, IgM < 1%, IgD < 1%
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.