



PRODUCT SPECIFICATION

CATALOG NUMBER: A37020

PRODUCT:	Human Immunoglobulin IgE Monoclonal Antibody
IMMUNOGEN:	Human IgE (Myeloma Serum)
FORM:	Protein A Purified immunoglobulin from mouse ascites in
PURIFICATION:	>90% by SDS-PAGE
BUFFER:	20 mM Phosphate, 150 mM Sodium Chloride, pH 7.2, pH 7.2
PRESERVATIVE:	0.05% NaN ₃
AFFINITY CONSTANT:	6.9 X 10 ¹⁰ L/m
SUBCLASS:	IgG ₁ , Kappa
STORAGE:	2-8 ^o C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.