



PRODUCT SPECIFICATION

CATALOG NUMBER: A37130

PRODUCT:	Human IgM Monoclonal Antibody
FORM:	Affinity purified immunoglobulin
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Mu-chain of human IgM
REACTIVITY:	Specific antibody titer is 1×10^6 by indirect ELISA. No cross-reactivity with human IgA, IgG, IgE
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.