



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A37140**

PRODUCT:	Human IgM Monoclonal Antibody
FORM:	Affinity purified immunoglobulin
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Mu-chain of human IgM.
REACTIVITY:	Specific antibody titer is $5 \times 10^{-7}$ by indirect ELISA. No cross-reactivity with human IgA, IgG, IgE
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.