



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A37150**

PRODUCT:	Human IgM Monoclonal Antibody
PRODUCT FORM:	Affinity purified mouse immunoglobulin
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity is tested by electrophoresis.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>2b</sub>
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Mu-chain of human IgM Fab region localized epitope
REACTIVITY:	Specific antibody titer is $5 \times 10^{-7}$ by indirect ELISA.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.