



PRODUCT SPECIFICATION

CATALOG NUMBER: A37160

PRODUCT:	Anti-human IgM Monoclonal Antibody
FORM:	Affinity purified immunoglobulin
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Mu-chain of human IgM, Fc region localized epitope
REACTIVITY:	Specific antibody titer is 5×10^{-7} by indirect ELISA.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.