



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A37190**

PRODUCT:	Human IgM Heavy Chain Monoclonal Antibody
FORM:	Protein A purified immunoglobulin
BUFFER:	PBS with 1 mg/ml bovine serum albumin.
PRESERVATIVE:	10mM sodium azide
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Human immunoglobulin M heavy chains.
APPLICATIONS:	Immunohistological staining of IgM expressing cells in tissue sections. Identification of surface IgM on B lymphocytes by indirect immunofluorescence. Identification of solid phase IgM by ELISA. Use as a solid phase "capture" antibody in sandwich ELISA techniques.