



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A37205**

PRODUCT:	IgM Monoclonal Antibody
IMMUNOGEN:	Human Serum IgM
LOT NUMBER:	Reference
PRODUCT FORM:	Purified IgG from Mouse Ascites
PURITY:	>95% by SDS-Page
BUFFER:	Phosphate Buffer Saline, pH 7.0
PRESERVATIVE:	0.02% Sodium Azide
SUBCLASS:	IgG <sub>1</sub>
CONCENTRATION:	2.17 mg/ml
STORAGE CONDITIONS:	Short term: 2-8°C Long term: Frozen
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.