



PRODUCT SPECIFICATION

CATALOG NUMBER: A37240

PRODUCT:	Human IgM Monoclonal Antibody
FORM:	DEAE purified immunoglobulin from mouse ascites
PURITY:	>90% (SPE gel showing one band only)
BUFFER:	10mM PBS, pH 7.2
PRESERVATIVE:	0.02% NaN ₃
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.