



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A37420**

PRODUCT:	IgG Fc Monoclonal Antibody
IMMUNOGEN:	Human IgG Fc
FORM:	Caprylic acid purified
PURITY:	>90% by SPE gel showing one band
BUFFER:	10mM PBS, pH 7.2
PRESERVATIVE:	0.02% NaN <sub>3</sub>
SUBCLASS:	IgG
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.