



## PRODUCT SPECIFICATION

|                      |  |
|----------------------|--|
| catalog number:      | <b>A38130</b>  |
| product description: | D-dimer Monoclonal Antibody  |
| lot number:          | [lot specific]   |
| source:              | Murine   |
| form:                | Liquid   |
| purification method: | Protein A chromatography   |
| purity:              | ≥90% by SDS-PAGE   |
| buffer:              | PBS, pH 7.4  |
| preservative:        | 0.02% NaN <sub>3</sub>   |
| isotype:             | IgG <sub>2a</sub>  |
| concentration:       | [lot specific] mg/mL   |
| pairing:             | This antibody can be used as the detection antibody in conjunction with A38140 (as the capture antibody).  |
| storage:             | 2-8°C short term (up to three months); -20°C long term (up to 36 months) avoiding multiple freeze thaw cycles.                                       |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known hazard/biohazard associated with this product. |