



## PRODUCT SPECIFICATION

catalog number:	<b>A38210</b>
product description:	Fibrin Degradation Product (FDP) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human d-dimer
source:	Murine
form:	Liquid
buffer:	0.015M sodium phosphate, 0.8% NaCl, pH 7.2
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This antibody recognizes human FDP and has cross-reactivity with human d-dimer.
pairing:	This antibody can be used in ELISA as the capture antibody with A38220 as the conjugate.
manufacture date:	[lot specific]
storage:	2-8°C short term (one year); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.