



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>A38220</b>   |
| product description: | Fibrin Degradation Product (FDP) Monoclonal Antibody, Purified  |
| lot number:          | [lot specific]  |
| immunogen:           | Human d-dimer   |
| source:              | Murine  |
| form:                | Liquid  |
| buffer:              | 0.015M potassium phosphate, 0.85% NaCl, pH 7.2  |
| preservative:        | 0.05% NaN <sub>3</sub>  |
| isotype:             | IgG <sub>1</sub>  |
| concentration:       | [lot specific] mg/mL  |
| specificity:         | This antibody recognizes human FDP and does not react with human d-dimer.   |
| pairing:             | This antibody can be used in ELISA as the conjugate antibody with A38210 as the capture.  |
| manufacture date:    | [lot specific]  |
| storage:             | 2-8°C short term (one year); -20°C long term avoiding multiple freeze thaw cycles.  |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |