



PRODUCT SPECIFICATION

CATALOG NUMBER: A40010

PRODUCT:	Prostatic Acid Phosphatase (PAP) Monoclonal Antibody
IMMUNOGEN:	PAP purified from seminal fluid
FORM:	Protein A purified immunoglobulin from mouse ascites
PURITY	>90% by SDS PAGE
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.05% Sodium Azide,
FILTER:	Filtered through 0.2um filter
AFFINITY CONSTANT:	1×10^{11}
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.