



PRODUCT SPECIFICATION

CATALOG NUMBER: A41050

PRODUCT:	Phenytoin Monoclonal Antibody
IMMUNOGEN:	Phenytoin-KLH
PURIFICATION:	Protein A column chromatography
BUFFER:	PBS, pH 7.0
PRESERVATIVE:	0.02% sodium azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.