



PRODUCT SPECIFICATION

CATALOG NUMBER: A42030

PRODUCT:	Phenobarbital Monoclonal Antibody
IMMUNOGEN:	Phenobarbital-BSA
PURIFICATION:	Protein A
PURITY:	>90% SDS PAGE
BUFFER:	10 mM Phosphate, 150mM NaCl, pH 7.4
PRESERVATIVE:	0.1% Sodium Azide
SUBCLASS:	IgG ₁ ,Kappa
STORAGE CONDITIONS:	2-8°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.