



PRODUCT SPECIFICATION

catalog number:	A42110
product description:	Phenobarbital Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Phenobarbital-BgG
source:	Murine
form:	Liquid
purification method:	Protein A
buffer:	0.05 M Phosphate, 0.85% sodium chloride, pH 7.2
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL
pairing:	This antibody forms a match pair with V55060 or V53070.
notes/applications:	ELISA
storage:	2-8°C short term (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.