



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A44020**

PRODUCT:	Human Prolactin Monoclonal Antibody
IMMUNOGEN:	Purified Human Prolactin
PURITY:	≥90% by protein A chromatography
BUFFER:	10mM Phosphate, 150 mM NaCl, pH 7.4
PRESERVATIVE:	.1% NaN <sub>3</sub>
AFFINITY CONSTANT:	9.3 X 10 <sup>9</sup>
SUBCLASS:	IgG1, Kappa
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.