



PRODUCT SPECIFICATION

CATALOG NUMBER: A44040

PRODUCT:	Human Prolactin Monoclonal Antibody
IMMUNOGEN:	Purified Human Prolactin
PURIFICATION:	Protein A, >90% pure by SDS-PAGE
BUFFER:	10mM Phosphate, 150 mM NaCl, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG _{2a} , Kappa
STORAGE CONDITIONS:	2-8°C short term, -20°C long term
COOPERATIVE PAIR:	Yes, with A44020
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.