



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A44090**

PRODUCT:	Human Prolactin Monoclonal Antibody
IMMUNOGEN:	Purified Human Prolactin
FORM:	DEAE column chromatography purified immunoglobulin
BUFFER:	0.015 M KPO <sub>4</sub> buffer, pH 7.40, containing 0.85% NaCl
PRESERVATIVE:	0.1% NaN <sub>3</sub>
STORAGE CONDITIONS:	2-8°C, do not freeze.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.