



PRODUCT SPECIFICATION

catalog number:	A44543
product description:	HIV-1 p24 Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for HIV-1 p24 antigen.
pairing:	This antibody can be used in sandwich ELISA as the conjugate antibody in conjunction with A44541 or A44542 (as the capture antibody).
storage:	2-8°C short term up to 2 years; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.