



## PRODUCT SPECIFICATION

catalog number:	<b>A44548</b>
product description:	HIV p24 Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Ion-exchange chromatography
purity:	>95%
buffer:	50 mM NaCl, 10 mM PBS, pH 8.0
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	HIV p24
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.