



PRODUCT SPECIFICATION

catalog number:	A44570
product description:	HIV-1 gp41 Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Ammonium sulfate precipitate
buffer:	PBS, pH 7.2
preservative:	0.01% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	HIV-1 gp41 antigen
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.