



## PRODUCT SPECIFICATION

catalog number:	<b>A44590</b>
product description:	HIV-2 gp36 Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein G
buffer:	PSB
preservative:	0.01% NaN <sub>3</sub>
isotype:	IgG <sub>2a</sub>
concentration:	[lot specific] mg/mL
specificity:	HIV gp36 antigen. Recognition epitope is aa 390-702 of HIV gp36 antigen.
application:	ELISA, Western Blot
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.