



PRODUCT SPECIFICATION

CATALOG NUMBER: A46040

PRODUCT:	Ca 15-3 Monoclonal Antibody
FORM:	Purified immunoglobulin from BALB/C mouse ascites
PURIFICATION:	Affinity chromatography
BUFFER:	15mM Phosphate buffer, 0.15M NaCl, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG _{1k}
STORAGE CONDITION:	2-8°C short term, -20°C for up to 12 months, avoid repeated freeze/thaws
SPECIFICITY:	Human Ca 15-3
PAIRING:	Can be used in Ca 15-3 sandwich ELISA as capture antibody in combination with A46050489P as conjugate antibody.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.