



PRODUCT SPECIFICATION

CATALOG NUMBER: A46074

PRODUCT:	CA 125 Monoclonal Antibody
FORM:	Purified immunoglobulin from BALB/C mouse ascites
PURIFICATION:	Affinity chromatography
BUFFER:	15mM Phosphate buffer, 0.15M NaCl, pH 7.2
PRESERVATIVE:	0.09% NaN ₃
SUBCLASS:	IgG _{1k}
STORAGE CONDITION:	2-8°C short term, -20°C for up to 12 months. Avoid repeated freeze/thaw cycles.
SPECIFICITY:	Human CA 125
PAIRING:	Can be used in CA 125 sandwich ELISA as the capture antibody in combination with A46075 or A46076 (conjugate).
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.