



PRODUCT SPECIFICATION

catalog number:	A46079
product description:	CA 125 Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Colorless liquid
purification method:	Affinity chromatography
buffer:	15mM phosphate buffer, 0.15M NaCl, pH 7.2
preservative:	0.09% NaN ₃
isotype:	IgG ₁ , kappa
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
specificity:	Human CA 125
application/notes:	This antibody can be used in CA125 sandwich ELISA as the detection antibody in conjunction with A46074 (as the capture antibody).
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.